


DOSSIER TECHNIQUE

SPECTROPHOTOMÈTRE BECKMAN DU 640

Comporte les documents suivants :

- | | |
|---|--------------|
| • Caractéristiques filtre interférentiel Ealing : | page 2 |
| • Caractéristiques filtres de densité Melles Griot: | page 3 |
| • Généralités sur les spectrophotomètres : | page 4 |
| • Configuration du Beckmann | page 5 |
| • Schéma optique DU640 : | page 5 |
| • Spécifications techniques DU640 : | page 6 |
| • Résumé notice : | page 7 et 8 |
| • Principe du "blanc" : | page 9 |
| • Modes d'acquisition des données : | page 10 à 15 |
| • Spectrophotomètre DU64 : | page 13 à 15 |




cadre 1 : Transmission du filtre interférentiel Ealing réf. 03FIM006.

MELLES GRIOT

Melles Griot BV
Aalsbergen 2 6942 SE Didam
Postbus 272 6900 AG Zevenaar The Netherlands

Tel +31 (0)316 333041
Fax +31 (0)316 528187
www.mellesgriot.com

ABN-AMRO Bank branch Zevenaar account no 43.80.96.177
Postbank branch Arnhem account no 31.12.96/
Chamber of Commerce Arnhem registration number 09072650



Description: 03FNG015 / OD=1,0

Spectrum Name: NC1001B.SP

Date Created: Tue May 04 09:39:17 2004

Instrument: Lambda 9 UV/VIS

Calibration Due Date: 25.09.2004

Serial Number # 101N1041605

Data Interval: 1.0000 nm

Scan Speed: 150 nm/min

Slit Width: 2.0000 nm

Smooth Bandwidth

 **Barloworld Scientific**



cadre 2 : Filtre de densité optique Melles Griot.

A. Généralités sur les spectrophotomètres

Un spectrophotomètre comprend 4 parties essentielles :


A.1. Source lumineuse

Elle est constituée par :

- Une lampe à décharge au deutérium utilisée dans le domaine de 190 à 400 nm avec un maximum d'émission à 652.1 nm (voir cadre 15).
- Une lampe à filament de tungstène pour la région allant de 350 à 800 nm (cadre 16).
- Parfois une lampe à décharge au xénon utilisée dans le domaine UV et visible. Ce type de lampe est très énergétique. Elle fonctionne sous forme de flash, juste au moment de faire une mesure.



cadre 3.



cadre 4.

A.2. Monochromateur

L'élément de base est un prisme, un réseau ou un filtre coloré. Le rôle du monochromateur est d'isoler le rayonnement sur lequel on fait la mesure. Il est composé principalement d'un système dispersif, d'une fente d'entrée et d'une fente de sortie (cadre 17).


A.3. Cuve

Elle contient soit l'échantillon, soit la référence. La longueur de la cuve est définie (1, 2, 4 ou 5 cm de trajet optique). Elle doit être transparente aux radiations d'étude. Par exemple en UV, les cuves sont en quartz, elles ne peuvent ni être en verre, ni en plastique.

A.4. Détecteur

A.4.1. Photodiode (semi-conducteur)


Lorsqu'un photon rencontre un semi-conducteur, il peut transférer un électron de la bande de valence (niveau énergétique bas) vers la bande de conduction (niveau énergétique haut) en créant une paire électron-trou. Le nombre de paires électron-trou est fonction de la quantité de lumière reçue par le semi-conducteur qui peut donc être utilisé en tant que détecteur optique.



cadre 5.

A.4.2. Barette CCD linéaire


L'emploi d'une barette de diodes permet une mesure simultanée sur toute l'étendue du spectre. Une barette CCD est un alignement de photodiodes de petites dimensions (14 x 14 μm) qui fonctionnent en intégrateur de lumière. La charge qui apparaît dans une photodiode est proportionnelle à l'exposition, c'est à dire au produit de l'éclairement par le temps de pose et elle dépend de la longueur d'onde. A la fin de la pose, le contenu des capteurs est transféré dans un registre analogique à décalage et une nouvelle pose commence. Ce registre transmet les données mémorisées en mode série, c'est à dire l'une après l'autre à un rythme fixé par l'électronique de commande de la barette CCD. Ces données apparaissent sous forme de grandeur électrique. Couplé à un ordinateur, le spectrophotomètre permet de tracer très rapidement des spectres d'absorption ou transmission. Le logiciel gère le temps de pose du capteur CCD.



cadre 6.


A.4.3. Photomultiplicateur

Une radiation (voir cadre 19) incidente arrache un électron de la cathode par effet photoélectrique. Cet électron est alors accéléré vers une seconde électrode appelée dynode portée à un potentiel supérieur. L'énergie de l'électron incident est suffisante pour arracher plusieurs autres électrons et ainsi de suite, d'où l'effet multiplicatif. Pour un électron arraché sur la cathode on peut récupérer une centaine d'électrons sur l'anode.




cadre 7.

B. Schéma optique



C. Configuration Beckmann

cadre 8 : Schéma optique du DU640.



cadre 9 : Configuration de la partie communication du Beckmann

D. Caractéristiques techniques du DU640

TECHNICAL SPECIFICATIONS

13.1 Performance Specifications

Scan Speeds	120, 240, 600, 1200, 2400 nm/min
Data Collection Rate	20 samplings/second
Response Time	0.05 second
Wavelength Range	190 to 1100 nm
Wavelength Accuracy at 656.1 nm	±0.2 nm
Full Range	±0.5 nm
Wavelength Repeatability at 656.1 nm	±0.1 nm
Full Range	±0.2 nm
Spectral Bandwidth (from 200 to 680 nm)	≤1.8 nm
Photometric Readout	-0.300 to 3.000 A or 0.0 to 200.0 %T
Photometric Accuracy (at 1 A with NIST 930D solid filter at 546 nm)	±0.005 A
RMS Noise	≤0.0002 A rms ¹ / ≤0.0005 A rms
(0A, average of 10 standard deviations of 10 readings at 0.05 sec intervals at 500 nm)	

cadre 10 : Spécifications techniques.

Stray Light (measured using NaI at 220 nm per ASTM E387-84)	<0.05%
Stability (0A, constant ambient conditions, measured for 1 hour at 340 nm)	<0.003 A
Baseline Flatness (from 200 to 900 nm, at 0 A)	±0.001 A rms ² / ±0.003 A rms
13.2 Physical and Environmental Specifications	
Width	69 cm (27 inches)
Height	58 cm (23 inches)
Depth	53 cm (21 inches)
Weight	33 kg (73 lbs)
Line Voltage	100/120V±10% or 220/240V±10%
Frequency	50/60 Hz
Power	200 watts typical
Ambient Temperature Operating Range	+15 to 40°C (59 to 104°F)
Humidity	<85% maximum relative humidity, not to exceed 32.5°C WBT

cadre 11 : Spécifications techniques.

E. Résumé notice DU640

Manual 517313

QUICK REFERENCE

DU® SERIES 600 SPECTROPHOTOMETER

REDIREAD

Use the RediRead™ mode to take absorbance readings on a sample at 280 nm.

Click on	Explanation
[UV OFF]	Turn on the UV source. Wait for source to light.
<<RediRead>>	Display the RediRead window.
486.0nm (or displayed wavelength value)	Display numeric keypad.
[2] [8] [0] [OK]	Change wavelength to 280 nm.
Insert blank or clear path (air blank).	
<ReadBlank>	Read and store blank at 280 nm.
Insert sample.	
<ReadSample>	Display absorbance reading at 280 nm.
<Print>	Print the window.
<Exit>	Remove the RediRead window.

REDISCAN

Use the RediScan™ mode to scan a sample.

Click on	Explanation
[VIS OFF]	Turn on the visible source.
[UV OFF]	Turn on the UV source.
<<RediScan>>	Display the RediScan window.
Insert blank or clear path (air blank).	
<ScanBlank>	Read and store the blank scan.
Insert sample.	
<ScanSample>	Read and display the sample scan.
<AutoScale> (absorbance limits)	Automatically scale the absorbance axis. (Change the absorbance limits, if desired.)
(wavelength limits)	Change the wavelength axis, if desired.
<Print>	Print the scan.
<Exit>	Remove the RediScan window.

WAVELENGTH SCAN

Use the Wavelength Scan mode to scan a sample in absorbance from 350 to 700 nm at 600 nm/min. Automatically scale the data and find the peaks. To set up the parameters:

- Verify that the Main window is displayed.
- Verify that the visible source is on.
- Click on "WAVELENGTH SCAN" to display the Wavelength Scan window.
- The following parameters should be listed near the top of the Wavelength Scan window. To change any parameter, click on the displayed value to display an input window. Input the desired value, then click on [OK] to remove the window.

Scan directory: VIEW
Start wl: 350 nm
End wl: 700 nm
Overlay scans: [No]
Autoprint: [No]
Autosave: [No]
Scans per sample: 1
Interval: 10.0 [sec]
Method name: A:\DEFAULT
Autosave name: A:\JSCANS
Sampling device: None
Scan speed: 1200 nm/min

- Click on <Function> to display the Function window. Click on to darken the boxes in front of "Scan" and "Peak Pick". Click on <Exit> to remove the Function window.

To analyze the samples:

Click on	Explanation
Insert blank or clear path (air blank).	
<<BLANK>>	Read and store the blank scan.
Insert sample.	
<ReadSamples>	Read and display the sample scan.
<AutoScale> (absorbance limits)	Automatically scale the data. Change the absorbance limits, if desired.
<Print>	Print the window.
<SaveClear>	Display the Save Clear window. To store the scan, click on the file name to display the alphanumeric keypad. Input the desired file name, then click on [OK], [OK] to store the data.
<Quit>	Display the Quit window. Click on [OK] to exit the mode.

FIXED WAVELENGTH

Use the Fixed Wavelength Scan to make absorbance readings at 230, 260, 280 and 320 nm. To set up the parameters:

- Verify that the Main window is displayed.
- Verify that the UV source is on.
- Click on "FIXED WAVELENGTH" to display the Fixed Wavelength window.
- Click on <Parameters> to display the Parameters window. To set up the desired parameters, the first four rows in the table should be the following:

Wavelength Factor	Units	Use
230	1.000	mg/ml
260	1.000	mg/ml
280	1.000	mg/ml
320	1.000	mg/ml

To change a parameter, click on the displayed value to display an input window. Input the desired value, then click on [OK] to remove the window. (Any information may appear in the other rows, except that [No] should be in the "Use" column.) When the desired information is displayed, click on <Exit> to remove the Parameters window.

- The following parameters should be listed near the top of the Fixed Wavelength window. To change a parameter, click on the displayed value to display an input window. Input the desired value, then click on [OK] to remove the window.

Results file: A:\WORK.RES
Read average time: 0.50 sec
Read mode: [Abs]
Method name: A:\DEFAULT
Sampling device: None

Click on	Explanation
Insert blank or clear path (air blank).	
<<BLANK>>	Read and store the blank readings.
Insert sample.	
<ReadSamples>	Read and display the sample readings.
→←	Slide the display horizontally to display data at all wavelengths.
<Print>	Print the window.
<SaveClear>	Display the Save Clear window. To store the data, click on the file name to display the alphanumeric keypad. Input the desired file name, then click on [OK], [OK] to store the data.
<Quit>	Display the Quit window. Click on [OK] to exit the mode.

cadre 12 : Résumé notice DU640.

TIME DRIVE	GRAPHIC MANIPULATIONS
Use the Kinetics/Time mode to monitor a reaction at 340 nm, every second for 1 minute. Plot the data with absorbance axis limits of 0.0 to 0.5. To set up the parameters: 1. Verify that the Main window is displayed. 2. Verify that the visible source is on. 3. Click on "KINETICS/TIME" to display the Plotting window. 4. The following parameters should be listed near the top of the Plotting window. To change any parameter, click on the displayed value for the parameter to display the numeric input window. Input the desired value, then click on [OK] to remove the window. Results file: A:\WORK RES Analytical wavelength: 340.0 nm Background wavelength: [No] 250.0 nm Sampling device: None Number of samples: 1 Sample assignment: 1 [S] Interval time: 1.000 [sec] Total time: 60.00 [sec] Method name: A:\DEFAULT Units: mg/ml Factor: 1.0000 Temperature: 25.0 C Curve offset: [None] Read average time: 0.50 sec To analyze the samples:	To change the axis limits: <u>Click on</u> <u>Explanation</u> <Autoscale> Automatically scale the ordinate axis. (any axis limit) Display numeric keypad. Click on desired axis value, then [OK] to remove the keypad and change the axis limit. <Zoom> To enlarge any portion of the graphics, click on two points on the graph to place crosses at the opposite corners of the area to be enlarged. When the second cross is clicked on, the graph is replotted. This can be repeated as often as desired. To return to the original plot, click on <ZoomOut>. To change y-axis label (Wavelength Scan, only): [Abs] Change the label to [Abs]. [%T] Change the label to [%T]. To display multiple functions on the y-axis (Wavelength Scan, only, without overlay): <Function> Display the Function window. Click on the desired function(s), then click on [OK] to remove the window and display the desired information. To find ordinate and abscissa values at any point on the graph: <Trace> Display the Trace window and place a vertical line on the graph. To move the vertical line on the graph, move the arrow to the point of interest on the graph and click on the center mouse button. To move the vertical line to either the right or left, click on the right or left mouse button, respectively. To remove the Trace window, move the arrow outside the graph and click on the center mouse button. To annotate the graph: <Annotate> Then position the arrow on the graph and click on the center mouse button to position a cross and display the alphanumeric keypad. Click on to input the desired information, then click on [OK].
Insert blank or clear path (air blank). <<BLANK>> Read the blank at 340 nm. <ReadSamples> Display Read Samples window. Insert sample. [START] Collect and display the data. When data collection is complete, the Read Samples window is removed. <Autoscale> Automatically scale the data. <Print> Print the window. <Quit> Display the Quit window. Click on [OK] to exit the mode.	BECKMAN

Beckman Instruments, Inc.
 2500 Harbor Boulevard • Fullerton, CA 92634-3100
 (714) 871-4848 • TWX: 010-922-1260 • Telex: 678413
 March 1992 © 1992 Beckman Instruments, Inc. Printed in U.S.A.

cadre 13 : Résumé notice DU640.

F. Réalisation d'un "blanc"

A blank is always required before data collection; any reading taken without a blank is invalid. A blank reading is taken when **< BLANK >** (located in the permanent menu bar on the bottom of the window) is clicked on.

NOTICE

In the RediRead Mode the blank command is **<ReadBlank>**. In the RediScan Mode, the blank command is **< ScanBlank>**.

When the instrument blanks, the following steps are performed:

1. The monochromator is moved to the proper wavelength. This is the specified wavelength for a single wavelength reading.
2. The proper detector gain value is selected automatically. This minimizes the noise level and maximizes photometric accuracy.
3. Dark current is measured and corrected. This compensation assures accurate readings at high absorbance.
4. In the Wavelength Scan mode, only, a background scan is made. The blank (or reference) is automatically scanned over the same range at the same speed that the samples will be scanned, so that the background correction is optimal.

This calibration assures repeatable readings every time the instrument is used.

In all modes, a blank solution should be in the sample compartment during the blank. It is suggested that the solvent used to prepare the samples be used for the blank. However, air (no sample) may be used. A new blank reading should be taken each time the solvent is changed.

NOTICE

Plastic cuvettes, glass (Pyrex) cuvettes, and some solvents have significant absorption in the UV region. Verify that they transmit UV light by scanning them versus air before using them in the UV region.

To re-zero the instrument at any time between samples, insert the same blank solution and click on **< BLANK >**.

The instrument stores the blank and uses it until either the sources are turned off or another blank reading is taken. For best results, the instrument should be blanked frequently, allowing the blank reading to be taken shortly before the sample measurement is taken. A new blank should be read if the instrument has not been used for an hour.

G. Modes d'acquisition des données

The DU Series 600 Spectrophotometer has five data collection modes: RediRead Mode, RediScan Mode, Fixed Wavelength, Wavelength Scan and Kinetics/Time.

G.1. RediRead Mode

The RediRead Window is used to take fixed wavelengths readings at one or more wavelengths quickly and easily. This window can be displayed whenever the instrument is not collecting data, regardless of the operating mode of the instrument. Data collected in this mode cannot be stored.

- 1) Click on **< RediRead >**, located in the permanent menu bar at the bottom of the display, to display the RediRead window, *figure 1*.
- 2) Sets the parameters :
 - a) Click on the wavelength value displayed and input the desired wavelength.
 - b) Click on **"Read avg time"** and input the desired read average time.
 - c) Verify that the desired reading mode is displayed, **[Abs]** or **[%T]**. Click on the mode to change it.

RediRead			HELP
ReadSample	ReadBlank	Print	Exit
500.0nm 0.1278 A			
Read avg time: 0.50		Read Mode: [Abs]	
Sample	Wavelength	Reading	
1	500.0nm	0.0109 A	
2	500.0nm	0.0481 A	
3	500.0nm	0.1091 A	
4	500.0nm	0.1278 A	
5			

- 3) Place a cuvette of solvent in the cell holder and click on **<ReadBlank>**. (If the instrument has previously been blanked at the selected wavelength using **< BLANK >**, it is not necessary to blank in the RediRead mode. **<ReadBlank>** in the RediRead mode does not affect the blank stored using **< BLANK >**.)
- 4) Place a cuvette of sample solution in the cell holder and click on **<ReadSample>**. The reading is displayed in the table on the window.
- 5) Repeat step 4 for all samples. The parameter input in step 2 can be changed at any time.

Readings from 11 samples are displayed on the window. When the sample 12 is read, the data is written over the data for sample 1.

- 6) To print the window, click on **<Print>**. Only the data that are displayed are printed.
- 7) To remove the RediRead window, click on **<Exit>**.

G.2. RediScan Mode

The RediScan window is used to make a wavelength scan at 1200 nm/min on a sample with minimum parameter setup. Data collected using this window cannot be stored ; the Wavelength Scan Mode must be used for data storage.

1. Click on **< RediScan >**, located in the permanent menu bar at the bottom of the display, to display the RediScan window, *figure 2*.
2. Verify that the proper ordinate label is displayed, **[Abs]** or **[%T]**. Click on the label to change it.
3. Verify that the desired wavelength limits are displayed. To change them, click on the displayed value and input the desired value. The sample will be scanned over the displayed wavelength range, only.
4. Place a cuvette of solvent in the cell holder and click **<ScanBlank>**. (If the instrument has previously been blanked in the Wavelength Scan mode at 1200 nm/min over the selected range, it is not necessary to blank in the RediScan mode).
5. Place a cuvette of sample solution in the cell holder and click on **<ScanSample>**. The scan data is displayed.
6. The following functions are available to reformat the data :
 - a) The data can be autoscaled by clicking on **<AutoScale>**.
 - b) Individual axis limit values can be changed by clicking on them and inputting the desired value.
7. To display the wavelength and ordinate readings at any point in the spectrum, click on **<Trace>**. Then move the mouse to the point of interest in the spectrum and click on the center mouse button to place a vertical line on the spectrum. The values at the place where the vertical line is placed are displayed in the lower right-hand

side of the window. To move the vertical line to either the right or left, click on the right or left mouse button, respectively.

8. To annotate the data, click on **<Annotate>**. Then click on the graph to position a cross and input information from the alphanumeric keypad or keyboard. Up to four annotations can be placed on the graph. The annotations are printed with the window, but are not stored with the data.
9. To print the wavelength scan in the window, click on **<Print>**.
10. Repeat steps 5 to 8 for all the samples.
11. To remove the RediScan window, click on **<Exit>**.

G.3. Fixed Wavelength

The Fixed Wavelength mode is used to collect data from a series of samples at up to 12 wavelengths. The data can be multiplied by user-input factor(s) to calculate a result at each wavelength. Any of the sampling devices can be used to simplify sample handling. Data can be stored for later recall.

To select the analysis parameters :

1. With the Main window displayed, click on **"FIXED WAVELENGTH"** to display the Fixed Wavelength window, *figure 3*.
2. Click on **<Parametres>** to display the Parameters window, *figure 4*.
 - a) Listed in the Parameters window are 12 wavelength values, with a factor and units that correspond to each wavelength. To change any of these values, click on the displayed value to display a keypad. Input the desired value on the keypad, then click on **[OK]** to accept the input and remove the keypad.
 - b) The fourth column in the Parameters window is the "Use" column. Each wavelength that is to be used in the analysis must have a "Yes" displayed. If a "No" is displayed for a desired wavelength, click on the "No" to display a "Yes".
 - c) When all the desired values are displayed, click on the **<Exit>** to remove the Parameters window. The input values are immediately displayed on the Fixed Wavelength window.
3. Readings can be taken either absorbance or transmittance. The selection is displayed following **"Read mode"** in the parameter listing near the top of the window. To change the read mode, click on the displayed option.




figure 2 : RediScan Window.

Fixed Wavelength							HELP
ReadSamples	Method	Parameters	SaveClear	Print	Quit		
Results file: A:\FIXED1		Method name: A:\FIXED		Sampling device: None			↑ + ↓
Read average time: 0.50		Read mode: [Abs]					
Sample ID	λ	Factor	λ	Factor	λ	Factor	
	350.0	56.00	440.0	230.0	520.0	6.500	
	Abs	Result	Abs	Result	Abs	Result	
		mg/ml		mg/ml		mg/ml	
1	0.2790	15.6221	0.1535	35.2955	0.3152	2.0490	
2	0.3647	20.4213	0.0971	22.3246	0.3784	2.4596	
43F	0.6747	37.7840	0.2244	51.6912	1.0832	7.0407	
43T	0.6413	35.9108	0.2421	55.6869	0.6864	4.4613	
46J	1.0447	58.5056	0.3162	72.7205	1.4303	9.2969	
48K	0.9504	53.2240	0.3767	86.6482	1.3800	8.9698	
7							

figure 3 : Fixed Wavelength Window.

Fixed Wavelength: Parameters				Exit
ClearAll	Print			
Wavelength	Factor	Units	Use	
350.0	56.00	mg/ml	[Yes]	
440.0	230.0	mg/ml	[Yes]	
520.0	6.500	mg/ml	[Yes]	
200.0	1.000	mg/ml	[No]	
400.0	1.000	mg/ml	[No]	
250.0	1.000	mg/ml	[No]	
300.0	1.000	mg/ml	[No]	
550.0	1.000	mg/ml	[No]	
600.0	1.000	mg/ml	[No]	
650.0	1.000	mg/ml	[No]	
700.0	1.000	mg/ml	[No]	
750.0	1.000	mg/ml	[No]	

figure 4 : Parameters Window.

To take readings :

1. Place a cuvette of solvent in the instrument. Click on < **BLANK** >.
2. If desired, click on the next displayed sample number and input up to an 11-gigit alphanumeric sample identification. If a sample identification is not input, the instrument numbers the samples consecutively.
3. Place a cuvette of sample solution in the cell holder and click on < **ReadSamples** >.
4. Data from up to 3 wavelengths are displayed at one time. To display data at other selected wavelengths, click on the right and left arrows, located on the right-hand side of the analysis parameters.
5. Repeat steps 2 to 4 until all samples have been read.
6. To print the sample data, click on < **Print** >.
7. When the analysis is complete, click on < **Quit** >. To store the method and/or results, click on the displayed file name(s) and input the desired file name(s). Then click on [OK] to store the data and return to the Main window.

G.4. Wavelength Scan

The Wavelength Scan mode is used to collect, manipulate and store scan data.

To select the analysis parameters :

1. With the Main window displayed, click on "**WAVELENGTH SCAN**" to display the Wavelength Scan window, figure 5.
2. Twelve parameters are listed near the top of the window.
 - a) Locate the "**Start wl**" and "**End wl**" parameters. To change the values, click on the displayed value and input the desired value.
 - b) Verify that the following parameters are as follow :

Overlay scans : [No]
 Autoprint : [No]
 Autosave : [No]
 Scans per sample : 1
 Sampling device : None
 Scan speed : 1200 nm/min




figure 5 : Wavelength Scan Window.

If any of the parameters are different than those listed, click on the displayed value and input the listed value.

3. The ordinate label and limits are displayed on the graphic portion of the window. To change any of these values, click on the displayed value and input the desired value.

To take readings :

1. Place a cuvette of solvent in the cell holder. Click on < **BLANK** >.
2. Place a cuvette of sample solution in the cell holder and click on < **ReadSamples** >.
3. The following functions are available to reformat the data :

Autoscale : Automatically scales the ordinate axis.

Limit changes : The limits on either axis can be changed by clicking on the displayed value and inputting the desired value.


Zoom : The "Zoom" feature is used to expand any portion of the graph. Click on < **Zoom** >, then click two points on the graph to place crosses at the opposite corners of the area to be enlarged. When the second cross is clicked on, the graph is replotted. This can be repeated as often as desired. To return to the original plot, click on < **ZoomOut** >.

4. To smooth the data, click on "**Smoothing**" and select the desired number of points of calculation. If too many points are used, real peaks can be lost.


5. To display a derivative or log scan, peak pick, valley pick and/or point pick, click on **<Function>** to display the Function Selection window. To select a function, click on to darken the box preceding the selection. Then click on **<Edit>** to remove the Function selection window from the display. The data are replotted using the selected function(s).
6. To display the wavelength and ordinate readings at any point in the spectrum, click on **<Trace>**. Then move the mouse to the point of interest in the spectrum and click on the center mouse button to place a vertical line on the spectrum. The values at the place where the vertical line is placed are displayed in the lower right-hand side of the window. To move the vertical line to either the right or left, click on the right or left mouse button, respectively.
7. To annotate, click on **<Annotate>**. Then, click on the graph to position a cross and input information from the alphanumeric keypad or keyboard. The annotation is printed with the window, but is not stored with the data.
8. To print the sample data, click on **<Print>**.
9. To store data before scanning another example, click on **<SaveClear>**. Click on the displayed file name, input the desired file name, then click on **[OK]**. The data are stored and the graphic area is cleared.
10. To scan additional samples, repeat steps 2 to 9, above.
11. When all the samples have been scanned, click on **<Quit>**. To store the method and/or displayed scan, click on the displayed file name(s) and input the desired file name(s). Then click on **[OK]** to store the data and return to the Main window.



cadre 14 : Spectrophotomètre DU64.



cadre 15 : Spectrophotomètre DU64.



cadre 16 : Spectrophotomètre DU64.