DOSSIER TECHNIQUE

SPECTROPHOTOMÈTRE SHIMADZU UV-1280

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Spectral bandwidth	5 nm
Wavelength range	190.0 to 1100.0 nm
Wavelength display	0.1 nm Increments
Wavelength setting	0.1 nm Increments (1 nm Increments when setting scanning range)
Wavelength accuracy	±1.0 nm
Wavelength repeatability	±0.3 nm
Wavelength slew rate	Approx. 6,000 nm/min Wavelength scanning speed: Approx. 9 nm/min to 1,600 nm/min
Light source switching	Select from the following three: •Automatic switching linked to the wavelength The switching wavelength can be selected between 295 nm and 364 nm, in 1 nm increments. The recommended wavelength is 340 nm. • Use the halogen (W1) lamp only, with no switching • Use the deuterium (D2) lamp only, with no switching
Stray light	0.05 % or less (220.0 nm NaI, 340.0 nm NaNO ₂)
Photometric system	Monitor double beam optics
Photometric range	Absorbance: -0.3 to 3.0 Abs Transmittance: 0.0 % to 200 %
Recording range	Absorbance: -4.0 to 4.0 Abs Transmittance: -400 % to 400 %
Photometric accuracy	±0.005 Abs at 1.0 Abs ±0.003 Abs at 0.5 Abs Using NIST 930D filter

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Photometric repeatability	±0.002 Abs at 1.0 Abs
Baseline stability	±0.001 Abs/h or less (700 nm, two hours after the power is turned ON)
Baseline flatness	±0.010 Abs or less (one hour after the power is turned ON, at 1,100 nm to 200 nm)
Noise level	P-P 0.002 Abs or less, RMS 0.0005 Abs or less
Baseline correction	Automatic correction via computer memory
Light source	20 W halogen lamp (2,000 hour operating life) Deuterium lamp (socket type) Built-in light source auto position adjustment
Monochromator	Uses an aberration correcting concave holographic grating
Detector	Silicon photodiode
Sample compartment	Internal dimensions: 110.0 (W) × 230.0 (D) × 105.0 (H) mm (Depth of one part is 155.0 mm.) Attachment method: two fastening screws
Display	6-Inch LCD (320×240 pixels) With LED illumination With contrast adjustment function
Output device	USB memory (optional) Data files saved in CSV format or UV-1280 dedicated format.
Power requirements	100 to 240 V, 50/60 Hz, 140 VA
Dimensions	W416 × D379 × H274 mm
Weight	10 kg
Environmental requirements	Temperature: 15 °C to 35 °C Humidity: 30 % to 80 % Humidity of 70 % or less at temperatures of 30 °C or higher

D. Principe du "blanc"

A blank is always required before data collection; any reading taken without a blank is invalid. A blank reading is taken when **< <BLANK> >** (located in the permanent menu bar on the bottom of the window) is clicked on.

NOTICE

In the RediRead Mode the blank command is **<ReadBlank>**. In the RediScan Mode, the blank command is **< ScanBlank>**.

When the instrument blanks, the following steps are performed:

1. The monochromator is moved to the proper wavelength. This is the specified wavelength for a single wavelength reading.
2. The proper detector gain value is selected automatically. This minimizes the noise level and maximizes photometric accuracy.
3. Dark current is measured and corrected. This compensation assures accurate readings at high absorbance.
4. In the Wavelength Scan mode, only, a background scan is made. The blank (or reference) is automatically scanned over the same range at the same speed that the samples will be scanned, so that the background correction is optimal.

This calibration assures repeatable readings every time the instrument is used.

In all modes, a blank solution should be in the sample compartment during the blank. It is suggested that the solvent used to prepare the samples be used for the blank. However, air (no sample) may be used. A new blank reading should be taken each time the solvent is changed.

NOTICE

Plastic cuvettes, glass (Pyrex) cuvettes, and some solvents have significant absorption in the UV region. Verify that they transmit UV light by scanning them versus air before using them in the UV region.

To re-zero the instrument at any time between samples, insert the same blank solution and click on **< <BLANK> >**.

The instrument stores the blank and uses it until either the sources are turned off or another blank reading is taken. For best results, the instrument should be blanked frequently, allowing the blank reading to be taken shortly before the sample measurement is taken. A new blank should be read if the instrument has not been used for an hour.

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