# DOSSIER TECHNIQUE

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## 6. DOCUMENTATION BRUKER

### 6.1. Advantages of FTIR spectroscopy


1) The sampling interval of the interferogram,  $\Delta x$ , is the distance between zero-crossings of the HeNe laser interferogram, and is therefore precisely determined by the laser wavelength. Since the point spacing in the resulting spectrum,  $\Delta \nu$ , is inversely proportional to  $\Delta x$ , FT-IR spectrometers have an intrinsically highly precise wavenumber scale (typically a few hundredths of a wavenumber). This advantage of FT spectrometers is known as **CONNES' advantage**.

2) The **JAQUINOT advantage** arises from the fact that the circular apertures used in FTIR spectrometers has a larger area than the slits used in grating spectrometers, thus enabling higher throughput of radiation.


3) In grating spectrometers the spectrum  $S(\nu)$  is measured directly by recording the intensity at successive, narrow, wavelength ranges. In FT-IR spectrometers all wavelengths from the IR source impinge simultaneously on the detector. This leads to the multiplex, or **FELGETT'S, advantage**.

The combination of the Jaquinot and Fellgett advantages means that the signal-to-noise ratio of an FT spectrometer can be more than 10 times that of a dispersive spectrometer.

#### IR spectrometer principle


Dispersive IR spectrometer

FT-IR spectrometer

### 6.2. The working principle of an FT-IR spectrometer

Infrared light emitted from a **source** (e.g. a SiC glower) is directed into an **interferometer**, which modulates the light. After the interferometer the light passes through the **sample compartment** (and also the sample) and is then focused onto the **detector**. The signal measured by the detector is called the **interferogram**.


General FT-IR spectrometer layout

### 6.3. Transmission spectrum


To calculate a transmission spectrum the following steps need to be performed :

- An interferogram measured without any sample in the optical path is Fourier transformed. This results in the single-channel reference spectrum  $R(\nu)$ .
- A second interferogram, measured with the sample in the optical path, is Fourier transformed. This results in the single-channel sample spectrum  $S(\nu)$ .  $S(\nu)$  looks similar to the reference spectrum, but shows less intensity at those wavenumbers where the sample absorbs radiation.

### 6.4. Michelson Interferometer

The **interferometer** is the heart of an FT-IR spectrometer. The collimated light from the infrared source impinges on a **beamsplitter**, which ideally transmits 50% of the light and reflects the remaining part. Having traveled the distance  $L$  the reflected light is hits a fixed mirror  $M1$ , where it is reflected and hits the beamsplitter again after a total path length of  $2L$ . The transmitted part of the beam is directed to a **movable mirror**  $M2$ . As this mirror moves back and forth around  $L$  by a distance  $Dx$ , the total path length is  $2(L + Dx)$ . The light returning from the two mirrors is recombined at the beamsplitter, with the two beams having a difference in path length of  $2Dx$ . The beams are spatially coherent and **interfere** with each other when recombined

### 6.5. Separation of spectral ranges

The mid-infrared, or MIR, is the spectral range from 4,000 to 400 cm<sup>-1</sup> wavenumbers. In this range fundamental vibrations are typically excited. In contrast, in the ‘near-infrared’, or NIR, spectral range, which covers the range from 12,500 to 4,000 cm<sup>-1</sup> wavenumbers, overtones and combination vibrations are excited. The far infrared’, or FIR, spectral range is between 400 and about 5 cm<sup>-1</sup> wavenumbers. This range covers the vibrational frequencies of both backbone vibrations of large molecules, as well as fundamental vibrations of molecules that include heavy atoms (e.g. inorganic or organometallic compounds).

The electromagnetic spectrum

Visible light and infrared light are two types of electromagnetic radiation, but with different wavelengths, or frequencies. In general, electromagnetic radiation is defined by the wavelength  $\lambda$  or the linear frequency  $\nu$ . The **wavelength** is the distance between two maxima on a sinusoidal wave.

The **frequency** is the number of wavelengths per unit time. Since all electromagnetic waves travel at the speed of light, the frequency corresponding to a given wavelength can be calculated as:

$$\nu = c/\lambda$$

According to the Plank’s Radiation Law, the frequency of electromagnetic radiation is proportional to its **energy**.

$$E = h \cdot \nu$$

$$\text{En consequence : } E = 1,24/\lambda \text{ (with } \lambda \text{ in } \mu\text{m} \text{ and } E \text{ in eV)}$$

In infrared spectroscopy **wavenumber** is used to describe the electromagnetic radiation. Wavenumber is the number of wavelengths per unit distance. For a wavelength  $\lambda$  in microns, the wavenumber, in cm<sup>-1</sup>, is given by

$$= 10000/\lambda$$

### 6.6. Fourier Transformation

Since the development of the first spectrophotometers in the beginning of 20th century a rapid technological development has taken place. The first-generation spectrometers were all dispersive. Initially, the dispersive elements were prisms, and later on they changed over to gratings. In the mid 1960s IR spectroscopy witnessed a revival due to the advent of spectrometers that utilized the Fourier transform (FT-IR). These second-generation spectrometers, with an integrated Michelson interferometer, provided some significant advantages compared to dispersive spectrometers. Today, almost every spectrometer used in mid-infrared spectroscopy is of the Fourier transform type. This is the reason why only FT-IR technology will be described in the following. Bruker Optics has specialized in the field of FT-IR spectroscopy since 1974, and is one of the leading manufacturers of FT-IR, FT-NIR and FT-Raman spectrometers throughout the world. The spectrometers are developed for analytical chemistry, life science, process, and many other fields.

Infrared spectroscopy is an extremely efficient analytical method due to modest operating expenditure. The analytical results are provided within a short period of time without the need of extensive sample preparation. In particular, infrared spectroscopy provides data which can be evaluated by quantity as well as by quality. The following will describe the qualitative and quantitative evaluation of acquired spectra.

- Qualitative evaluation of spectra
  1. Identify an unknown substance
  2. Check the identification of a known substance
- Quantitative evaluation of spectra

### 6.7. Transmission spectrum

To calculate the transmission spectrum the following steps need to be performed:

- an interferogram measured without any sample in the optical path is Fourier transformed. This results in the so-called single-channel reference spectrum  $S(\nu)$ .
- A second interferogram, measured with the sample in the optical path, is Fourier transformed. This results in the single-channel sample spectrum  $S(\nu)$ .  $S(\nu)$  looks similar to the reference spectrum, but shows less intensity at those wavenumbers where the sample absorbs radiation.

The final transmission spectrum  $T(\nu)$  is obtained by dividing the sample spectrum by the reference spectrum:

$$T(\nu) = S(\nu)/R(\nu)$$

### 6.8. Quantitative evaluation of spectra

The basic principle for quantitative evaluation in optical spectroscopy as well as in IR spectroscopy is the **Bouguer-Lambert-Beer Law** which had already been defined in 1852. Quantitative determinations by means of IR spectroscopy are preferably performed in solution. **Transmission T** of a sample is defined as:

$$T = I / I_0$$

$I_0$  is the intensity of the incident light beam,  $I$  is the intensity of the light beam leaving the sample. The percentage transmission (%T) is  $100 \cdot T$ . When traversing the measurement cell, the light intensity decreases exponentially:

$$I = I_0 \cdot \exp(-2.303\epsilon \cdot c \cdot b)$$

Where  $\epsilon$  is the **molar absorption coefficient** (in  $L \cdot mol^{-1} \cdot cm^{-1}$ ),  $c$  is the sample **concentration** (in  $mol \cdot L^{-1}$ ) and  $b$  the **thickness** of the measurement cell (in cm). The absorption coefficient  $\epsilon$  is a value which depends on either the wavelength or the wavenumber, which is typical for the compound analyzed. From the equation above, it follows that:

$$\log(I / I_0) = -\epsilon \cdot c \cdot b, \text{ or:}$$

$$A = \log(I_0 / I) = \epsilon \cdot c \cdot b$$

where  $A$  is the **absorbance**. Because of the **Bouguer-Lambert-Beer Law**, the relationship between **absorbance** and **concentration** of the absorbing substance is a **linear** function.